



Particle Size Analyzer (PSA)

Particle Size Analyzer covers a broad field of applications in practically any industry. Whenever quality control of powder materials or emulsions is needed, PSA offers a quick and reliable solution to monitor the particle size distribution.

Compact 2-in-1 design to measure both wet and dry samples

Multi-laser technology for a broad measurement range

- ❖ Dry Samples: 0.1 μm to 2500 μm
- ❖ Wet Samples: 0.04 μm to 2500 μm

For enquiries, please contact the staff-in-charge.

Name: Mr. Tan Han Kiang

Phone: 67905350

Email: chktan@ntu.edu.sg

Usage Rate: [Key Equipment Charges For ELC.xlsx](#)